



EXCEPTIONAL PRODUCTIVITY IN A COMPACT DESIGN

A powerful benchtop system that
maximizes valuable laboratory space

Access 2 Immunoassay System



Request a quote. 

MEDTEK is a distributor of Beckman Coulter.

200 J. Elizalde St. BF Homes, Parañaque City

Landline: +632 8806 9267, +632 8806 9943

Email : info@medtek.com.ph

Website: <https://medtek.com.ph>



Move healthcare forward.



ACCESS IMMUNOASSAY MENU

Adrenal/Pituitary

- > Cortisol (Serum/Plasma/Urine)

Anemia

- > Active B12[†]
- > EPO
- > Ferritin
- > Folate
- > Folate (RBC)
- > Intrinsic Factor Ab
- > sTfR
- > Vitamin B12

Blood Virus

- > HAV Ab
- > HAV IgM
- > HBc Ab
- > HBc IgM
- > HBs Ab
- > HBs Ag
- > HBs Ag Confirmatory
- > HCV Ab V3[‡]
- > HIV Combo[‡]

Bone Metabolism

- > Intact PTH
(Routine or Intraoperative)
- > Ostease
- > Ultrasensitive hGH
- > 25(OH) Vitamin D Total

Cardiovascular

- > CK-MB
- > Digoxin
- > hsTnI
- > Myoglobin
- > Triage BNP[§]

Diabetes

- > Ultrasensitive Insulin

Infectious Disease

- > CMV IgG
- > CMV IgM
- > Rubella IgG
- > Rubella IgM
- > Toxo IgG
- > Toxo IgM II

Inflammation

- > IL-6 (RUO in Japan)

Reproductive

- > AFP (ONTD)
- > AMH
- > DHEA-S
- > hFSH
- > hLH
- > Inhibin A
- > PAPP-A
- > Progesterone
- > Prolactin


- > Sensitive Estradiol
- > SHBG
- > Testosterone
- > Total β hCG (5th IS)
- > Unconjugated Estrinol

Thyroid

- > T3 (Free and Total)
- > T4 (Free and Total)
- > Thyroglobulin
- > Thyroglobulin Ab II
- > Thyroid Uptake
- > TPO Ab
- > TSH (3rd IS)

Tumor Markers

- > AFP
- > CEA
- > BR Monitor (CA 15-3 Antigen)
- > GI Monitor (CA 19-9 Antigen)
- > OV Monitor (CA 125 Antigen)
- > Prostate Health
Hybritech PSA
Hybritech free PSA
Hybritech p2PSA
Prostate Health Index (*phi*)^{||}



Access Immunoassay Menu

Beckman Coulter offers an extensive immunoassay diagnostic testing menu to fulfill patient-care needs.

- > More than 50 diagnostic immunoassays available
- > Multiple assays available for different diagnostic specialties
- > Continuous menu expansion with new and novel assays

Easy-to-use, easy-to-store reagent kits

The Access 2 uses liquid, ready-to-use reagent packs that contain everything needed to run an assay in one convenient container.

- > Slimline design facilitates easy pack loading
- > Self-sealing elastomer closure
- > No caps, preventing potential probe crashes
- > Easy-to-store kit boxes maximize refrigerator space
- > Typical shelf-life of 12 months

* TSH (3rd IS) is available in 100-test pack.

[†] In development.

[‡] Available through a partnership with Bio-Rad Laboratories, Inc.

[§] Available exclusively from Quidel or its authorized distributors for use on Beckman Coulter immunoassay systems. Quidel and Triage are trademarks of the Quidel group of companies.

^{||} *phi* is a combination of Hybritech PSA, free PSA and p2PSA.

Brochure is not for distribution in the U.S.

© 2018 Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries. For Beckman Coulter's worldwide office locations and phone numbers, please visit www.beckmancoulter.com/contact

BR-45565-v3 WF#BR-64595



Move healthcare forward.