

ARIA RAPID TESTS

PRODUCT CATALOGUE



DENGUE IgG/IgM COMBO RAPID TEST

Intended Use: Qualitative detection and differentiation of IgG and IgM Antibodies against Dengue Virus in human serum, plasma or whole blood

Reading Time: 25 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



DENGUE Ag RAPID TEST

Intended Use: Qualitative detection of Dengue Antigen (Ag) in serum, plasma or whole blood

Reading Time: 20 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



DENGUE DUO (Ag-IgG & IgM) RAPID TEST

Intended Use: Qualitative detection and differentiation of IgG and IgM Antibodies and Dengue Antigen (Ag) against Dengue Virus in human serum, plasma or whole blood

Reading Time: 20 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



CHIKUNGUNYA IgM COMBO RAPID TEST

Intended Use: Qualitative detection of IgM Antibodies against Chikungunya Virus in human serum, plasma or whole blood

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



HAV IgM RAPID TEST

Intended Use: Qualitative detection of IgM Antibodies against Hepatitis A Virus (HAV) in human serum or plasma

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



TYPHOID IgG/IgM COMBO RAPID TEST

Intended Use: Qualitative detection of IgG and IgM Antibodies against Salmonella typhi and Paratyphi in human serum, plasma or whole blood

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



H. PYLORI Ab COMBO RAPID TEST

Intended Use: Qualitative detection of IgG, IgM and IgA Antibodies against Helicobacter Pylori in human serum, plasma or whole blood

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



H. PYLORI Ag COMBO RAPID TEST

Intended Use: Qualitative detection of antigen against Helicobacter Pylori in human stool

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 25 tests
--------------------------	-----------------------	--------------------------



LEPTOSPIRA IgG/IgM COMBO RAPID TEST

Intended Use: Qualitative detection of IgG and IgM Antibodies against Leptospira interrogans in human serum, plasma or whole blood

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 30 tests
--------------------------	-----------------------	--------------------------



FECAL OCCULT BLOOD (FOB) RAPID TEST

Intended Use: Qualitative detection of fecal occult blood in human stool

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 25 tests
--------------------------	-----------------------	--------------------------



STREP A (DIP STRIP) RAPID TEST

Intended Use: Qualitative detection of Streptococcus Group A Antigen using throat specimen

Reading Time: 15 minutes	Shelf Life: 18 months	Packaging Size: 25 tests
--------------------------	-----------------------	--------------------------