



Coag2D

THE MOST VERSATILE 2 CHANNEL OPTICAL COAGULOMETER AT YOUR HAND

small compact design
supports D-dimer and AT-III test
ready for coagulation, chromogenic and turbidometric tests



MEDTEK

Request a quote. 

MEDTEK is a distributor of Diagon products.
200 J. Elizalde St. BF Homes, Parañaque City
Landline: +632 8806 9267, +632 8806 9943
Email : info@medtek.com.ph
Website: <https://medtek.com.ph>



Coag 2D

COAG 2D COAGULOMETER SPECIFICATIONS

Key features:

- 2 independent measuring channels
- 2 reagent positions (2 mL, 5mL)
- 4 incubation positions
- Coagulation screening test
- Quantitative D-dimer measurement
- AT-III activity measurement
- QC feature
- Supports application of Master Calibration Curve
- Maintenance free
- Up to 1000 results storage
- External data archiving possibility
- Connection to Coag Data data management software

System elements:

- 2 tempered reagent positions with magnetic stirring for heating and mixing
- Colour LCD touchscreen
- Innovative user interface
- LIS interface connection
- Small, compact footprint, easy to operate
- 405 nm, 570 nm and 640 nm wavelength detection
- Optional BCR and printer



EUROPE

Diagon Ltd - Hungary
1047 Budapest, Baross u. 48-52.
phone.: (+36 1) 369 6500
fax: (+36 1) 369 6301
e-mail: diagon@diagon.com

Diagon GmbH - Austria
GZO - Zukunftsweg 1,
Siegendorf, Burgenland 7011
phone: (+43 2687) 42090
e-mail: diagonmbh@diagon.com
www.diagon.com

SOUTH AMERICA

Diagon Ltda - Brasil
Av. Joaquim José Diniz, 787
CEP 31910-520
Bairro Fernão Dias
Belo Horizonte/MG
phone: (+55 31) 3488 6777
fax: (+55 31) 3488 6778
e-mail: diagon@diagon.com.br
www.diagon.com.br

ASIA

PT Inti Diagonama Selaras - Indonesia
Perkantoran Plaza Pasifik A3/59 -
Jl. Raya Boulevard Barat
Kelapa Gading Permal Jakarta 14250
phone.: (+62 21) 458 40649, 458 40650
fax: (+62 21) 458 45024
e-mail: dsu@centrin.net.id