



**MEDTEK**<sup>TM</sup>  
Confidence in Medical Diagnostics

**QA** QUALITY ASSURANCE **ISO 9001:2015**



# LAURA M

AUTOMATED URINE ANALYZER FOR  
PRECISE PATIENT DIAGNOSIS

- Up to 600 strips per hour
- High memory capacity
- Built-in waste container for used test strips

# LAURASMART

THE SMART SOLUTION FOR ACCURATE  
URINE ANALYSIS

- Automatic recognition of PHAN LAURA test strips
- High throughput (up to 240 tests per hour)
- User-friendly operation menu



**Erba**<sup>®</sup>  
Mannheim



# LAURA DEKAPHAN

URINE TEST STRIPS FOR  
COMPLETE URINE PANEL TESTING

10 parameters: Leukocytes, Erythrocytes, Urobilinogen, Bilirubin, Glucose, Protein, pH, Specific Gravity, Nitrite, Ketones

# LAURA MICROALBUPHAN

URINE TEST STRIPS FOR  
DETERMINATION OF MICROALBUMINURIA

- Quantitative determination of microalbumin
- Quantitative determination of creatinine
- Automatic calculation of albumin/creatinine ratio



# LAURA TETRAPHAN

URINE TEST STRIPS FOR  
ROUTINE URINE ANALYSIS

- Quantitative determination of glucose
- Quantitative determination of protein
- 4 parameters: Glucose, Protein, pH, Specific Gravity

## URINE STRIPS



### WITH ASCORBIC ACID PROTECTION

The PHAN LAURA urine strips are equipped with Ascorbic Acid Protection which prevents excessive presence of Ascorbic Acid in urine — ensuring accurate & reliable results.